

CS463/516

Medical Imaging


Lecture 4

Warning – theoretical lecture, mostly physics.


We will get to more traditional image processing topics soon 😊

Magnetic Resonance Imaging (MRI)




- Protons and neutrons in nucleus of atom possess angular momentum called *spin*
- Can observe spin in atoms with odd number of nuclei (such as hydrogen)
 - In atoms with even number of nuclei, the spins cancel each other out
- MRI machine produces static magnetic field B_0 , field strength measured in *Tesla*
- B_0 causes spins to align **parallel** or **antiparallel** to magnetic field direction
- Spin frequency ω depends on:
 - 1) strength of B_0
 - Measured in Tesla, typically 3 Tesla for clinical scanner
 - 2) atom-specific *gyromagnetic constant* γ
 - Measured in *MegaHz/Tesla*
 - $\omega = \gamma B_0$. Also known as the *Larmor frequency*
 - Example: hydrogen gyromagnetic ratio = $42.58 \frac{\text{MHz}}{\text{Tesla}}$
 - $\Rightarrow 3 \text{ Tesla} * 42.58 \frac{\text{MHz}}{\text{Tesla}} = 127.74 \text{ MHz}$
 - We call this the *resonant frequency* of hydrogen at 3 Tesla



Free induction decay (FID)




- How to measure a signal from all these spinning protons?
 - A spinning proton is a moving charge, which creates a magnetic field (electromagnetic induction)
 - However, spins are aligned with B_0 which is orders of magnitude stronger (signal from spins is obscured)
 - We therefore use a *radiofrequency (RF) pulse* perpendicular to B_0 to tilt spins away from B_0 (a)
 - RF pulse is an electromagnetic wave we construct to match hydrogen resonant frequency ω
 - Switch on RF pulse which deposits energy, tilting all spins by angle $0^\circ < \alpha < 180^\circ$
 - Switch off the RF pulse, spins will slowly return to original direction while still precessing (b)
 - c) Can now observe spins because they are out of alignment with B_0 . Observed signal is called *free induction decay*



Longitudinal and transverse magnetization (T1 and T2)


- After RF pulse tilts spins, they precess in transverse (xy) plane, gradually re-aligning with B_0 (a)
 - This process is called *relaxation*
- Therefore, we speak of two magnetization components (b)
 - 1) **longitudinal** magnetization or M_z
 - 2) **transverse** magnetization or M_{xy}
- Biological tissue can be characterized by two different relaxation mechanisms:
 - **T1 relaxation** – the process by which net magnetization returns parallel to B_0
 - **T2 relaxation** – the process by which transverse component of magnetization decays




MRI image contrast

- What does it mean to say: ‘an image is T1-weighted’?
 - It means, roughly speaking, we are measuring in each voxel how long it takes for longitudinal magnetization (M_z) to recover
 - Tissue types where M_z recovers quickly will appear bright on a T1-weighted image – example: the *white matter* which is mostly fat (a)
- What does it mean to say: ‘an image is T2-weighted’?
 - It means, roughly speaking, we are measuring in each voxel how long it takes for transverse magnetization M_{xy} to decay
 - Tissue types where M_{xy} decays slowly (long T2) will appear bright on T2-weighted image – example: cerebrospinal fluid which is mostly water (b)


Should measure signal at this time for optimal T1 contrast



T1-weighted




T2-weighted



Closer look at T2 relaxation



- a) the fluids appear bright on a T2-weighted image because transverse magnetization (M_{xy}) decays more slowly in fluids than in solid tissue
 - M_{xy} decay also called 'spin-spin relaxation'
- Can take 3-4 seconds for M_{xy} to decay in water
- Why does M_{xy} decay more slowly in fluid?
- Reasons for M_{xy} decay are complex, but basically, in water, M_{xy} decay occurs more slowly because spin-spin interactions aren't as strong as in other tissues
- b) M_{xy} decays exponentially with time due to *dephasing* (more on this later)




Why are we spending so much time on MRI? Because most of the images you will see in this course will be MRI, so its useful to be able to interpret image intensities.

Closer look at T1 relaxation

- a) solid tissue (fat in this example) appear bright in a T1-weighted image because longitudinal magnetization (M_z) recovers more quickly in solid tissues
 - T1 relaxation also called 'spin-lattice relaxation'
- Why does M_z recover more quickly in solids than in fluids?
- M_z recovery is fastest when molecular tumbling rate is equal to Larmor frequency
- Water has wide range of tumbling rates, so M_z recovery is slower in water
- b) T1 relaxation is the process by which net magnetization returns to its initial maximum value M_0





Magnetic resonance imaging



- We saw basic MRI experiment is as follows:
 - **1)** place sample in strong magnetic field B_0 , causing hydrogen atoms (spins) to align with B_0
 - **2)** apply RF pulse to deposit energy into sample, tilting spins away from B_0 direction
 - If you want to understand *why* this works, take PHY361 Quantum Mechanics
 - **3)** spins precess back into alignment with B_0 , creating measurable signal (electromagnetic induction)
 - **4)** decide some optimal time to measure signal from this changing magnetic field
 - Time at which we measure is called 'echo time' or TE, we set TE to achieve maximum contrast
- Problem? This doesn't give us an **image**.
 - RF pulse tilts spins in entire sample (all spins have same resonant frequency due to B_0)
 - This basic experiment captures signal from *entire sample*. We want to know signal from each voxel!
 - We need a way to *spatially encode* which spins our RF pulse will tilt away from B_0
 - Solution? Magnetic field *gradients* (next slide)


Slice selection with magnetic field gradients

- Recall the Larmor frequency: $\omega = \gamma B_0 = 127.74$ MHz for hydrogen at 3 Tesla
 - Basically, all hydrogen atoms are spinning at 127.74 MHz, if placed in a 3 T magnetic field
- To *select a slice* for imaging, apply *magnetic field gradient* (a)
 - Changes slightly the Larmor frequency of spins (gradient adds to B_0)
 - Example (a) – spins slightly slower near feet, spins slightly faster near head
 - Protons now spin with frequency that depends on their location in Z
- Craft RF pulse to match Larmor frequency of slice we want to image
 - Typically, *sinc pulse* (b) is used because it's FFT is a rectangle (allows to select rectangular slice)
 - Applying sinc pulse while gradient is on allows to selectively measure signal from spins in a single slice of the body we are attempting to image




Frequency encoding

- Gradient applied and sinc pulse delivered. Now measuring spins from single slice.
- How to localize spins within a slice? Use *frequency encoding*
- Frequency encoding:
 - Turn off original (Z) gradient used for slice selection
 - Apply gradient in x-direction, recorded signal now contains a *range* of frequencies (a)
 - Frequencies w_0 to w_{max} correspond to x-values x_0 to x_{max} .
 - Do FFT of measured signal, magnitude of frequency $w_0 + (w_{max} - w_0) \cdot (x_k - x_0) / (x_{max} - x_0)$ now corresponds to response of all protons on line $x = x_k$
 - This produces projection of spin density along y onto the x axis
 - Complete set of projections can be obtained by adding projections at measurements with gradient fields along different directions.
 - Using this data, can reconstruct image using filtered backprojection
 - This was the original method. Now we use k-space imaging (next slide)




K-space imaging (frequency space imaging)

- K-space imaging: generate image in frequency space, then apply ifft to get image
- 1) apply gradient G_z in z-direction to isolate a slice
- 2) apply another linear gradient in y-direction, increasing spin frequency along y, spins dephase (this is called the *phase encoding gradient*)
- 3) apply frequency encoding gradient and acquire data, signal is transformed into frequency space
- Result is integrated spin densities along x weighted with cosine wave of frequency 1 for each line of constant y.
- In other words, we have produced real-valued part of next line in k-space of the Fourier transform in the k_x direction. Imaginary part is computed by weighting signal with corresponding sine wave.
- Repeat the measurement with phase encoding for remaining lines in k-space. Once k-space is filled, image is generated by transforming it back into spatial domain (ifft)



K-space (frequency space)

Each cell in the grid represents a 2d wave with specific spatial frequency and orientation



Signal Generation and Spatial Encoding



David M. Higgins PhD

2019-05-12 ISMRM Montreal

- Excellent talk by David Higgins. Watch this talk, come back and re-read slides. Repeat until you understand.
- <https://youtu.be/29iffcZPQeM>